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Edited by M. Yamada, 8/21/2007

Lectin No.	Lectin	Reported specificity
1	LTL	Fuca1-3(Galb1-4)GlcNAc, Fuca1-2Galb1-4GlcNAc
2	PSA	Fuca1-6GlcNAc, a-D-Glc, a-D-Man
3	LCA	Fuca1-6GlcNAc, a-D-Glc, a-D-Man
4	UEA-I	Fuca1-2Galb1-4GlcNAc
5	AOL	Fuca1-6GlcNAc (core fucose)
6	AAL	Fuca1-6GlcNAc, Fuca1-3(Galb1-4)GlcNAc
7	MAL	Siaa2-3Galb1-4GlcNAc
8	SNA	Siaa2-6Gal/GalNAc
9	SSA	Siaa2-6Gal/GalNAc
10	TJA-I	Siaa2-6Gal/GalNAc
11	PHAL	tri/tetra-antennary complex-type N-glycan
12	ECA	Galb1-4GlcNAc
13	RCA120	Galb1-4GlcNAc
14	PHAE	bi-antennary complex-type N-glycan with outer Gal and bisecting GlcNAc
15	DSA	(GlcNAcb1-4) <sub>n</sub> , Galb1-4GlcNAc
16	GSL-II	agalactosylated tri/tetra antennary glycans, GlcNAc
17	NPA	High-Mannose, Mana1-6Man
18	ConA	High-Mannose, Mana1-6(Mana1-3)Man
19	GNA	High-Mannose, Mana1-3Man
20	HHL	High-Mannose, Mana1-3Man, Mana1-6Man
21	ACG	Siaa2-3Galb1-4GlcNAc
22	TxLCI	Mana1-3(Mana1-6)Man, bi- and tri-antennary complex-type N-glycan, GalNAc
23	BPL	Galb1-3GalNAc, GalNAc
24	TJA-II	Fuca1-2Galb1-> or GalNAcb1-> groups at their nonreducing terminals
25	EEL	blood group B antigen, Gala1-3Gal
26	ABA	Galb1-3GalNAc
27	LEL	GlcNAc trimers/tetramers
28	STL	GlcNAc oligomers, oligosaccharide containing GlcNAc and MurNAc
29	UDA	GlcNAcb1-4GlcNAc, Mixture of Man5 to Man9
30	PWM	(GlcNAcb1-4) <sub>n</sub>
31	Jacalin	Galb1-3GalNAc, GalNAc
32	PNA	Galb1-3GalNAc
33	WFA	GalNAcb1-4GlcNAc, Galb1-3(-6)GalNAc
34	ACA	Galb1-3GalNAc
35	MPA	Galb1-3GalNAc, GalNAc
36	HPA	a-linked terminal GalNAc
37	VVA	a-linked terminal GalNAc, GalNAca1-3Gal
38	DBA	blood group A antigen, GalNAca1-3GalNAc
39	SBA	a- or b-linked terminal GalNAc, GalNAca1-3Gal
40	Calsepa	Mannose, Maltose
41	PTL-I	a-linked terminal GalNAc
42	MAH	Siaa2-3Galb1-3(Siaa2-6)GalNAc
43	WGA	chitin oligomers, Sia
44	GSL-I A4	a-linked GalNAc
45	GSL-I B4	a-linked Gal

Remarks) These data were collected from lectin vendors and reports found by internet searches.

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